

# Security with every Stitch. Savings with every Suture.

From the world's largest independent suture manufacturer, Surgical Specialties



Quality veterinary sutures,  
40 years in the making.

TruStitch sutures offer remarkable value,  
quality, and precision that's comparable  
to significantly more expensive brands.

## ✓ Superior Needle Performance

**German-Engineered:** Characterized by superior sharpness, penetration and strength

**Consistency & Repeatability:** Manufactured with automated, robotic technology







## ✓ Clear Field of View

**Silver or Midnight Black Needles:** Midnight Black low-glare suture needles offer enhanced visibility in the surgeon's field of view.

TRUSTITCH™  
Veterinary Sutures

# Available in the leading animal health combinations.

## Absorbable Suture

MATERIAL TYPE		IN VIVO STRENGTH RETENTION	MASS ABSORPTION
<b>MONODERM</b> (Polyglycolic acid and e-caprolactone)		42% to 76% at 7 days 36% to 52% at 14 days	Essentially complete by <u>90 days</u>
<b>PDO</b> (Polydiocanone)		80% to 90% at 14 days 60% to 82% at 28 days 47% to 79% at 42 days	Essentially complete between <u>180-220 days</u>
<b>POLYSYN</b> (Polyglycolic acid)		50% at 14 days 20% at 21 days	Essentially complete between <u>50-90 days</u>
<b>CHROMIC GUT</b>		Varies	Varies

### FEATURES:

Dyed suture

USP sizes



6-0 through 1

18" • 30" • 36"  
lengths

Race track packaging

12 per box

## Non-Absorbable Suture

MATERIAL TYPE		IN VIVO STRENGTH RETENTION	MASS ABSORPTION
<b>NYLON</b>		Permanent	Permanent
<b>POLYPROPYLENE</b>		Permanent	Permanent



TruStitch Quill Vet Ethicon™ Covidien™

### Precision Reverse Cutting Needles



3/8 Circle

P-3	P-3	P-13	13mm
PS-2	PS-2	P-12	19mm
PS-1	PS-1	P-14	24mm
PS	PS		26mm

### Reverse Cutting Needles



3/8 Circle

FS-2*	FS-2	FS-2	C-13	19mm
FS-1*	FS-1	FS-1	C-14	24mm
FSL	FSL	FSL	C-16	30mm
FSLX*	FSLX	FSLX	C-17	36mm

1/2 Circle

CP-2*	CP-2	CP-2	GS-10	26mm
CP-1*	CP-1	CP-1	GS-11	36mm
CP*	CP	CP	GS-12	40mm

TruStitch Quill Vet Ethicon™ Covidien™

### Taper Point Needles



1/2 Circle

RB-1*	RB-1	RB-1	CV-23	17mm
SH-1*	SH-1	SH-1	CV-25	22mm
SH*	SH	SH	V-20	26mm
CT-2*	CT-2	CT-2	GS-22	26mm
CT-1*	CT-1	CT-1	GS-21	36mm

### Straight Cutting Needles



KS	KS	SC-2	60mm
----	----	------	------



\*Available in low-glare  
TruStitch Black Needles

Visit us at [surgicalspecialties.com](https://www.surgicalspecialties.com) for complete listing of needle and suture combinations.

**TRUSTITCH**  
Veterinary Sutures