

Choose Excellence,
Expect Value

SharpPoint  **PLUS**
Surgical Sutures

Suture Catalog



**SURGICAL
SPECIALTIES**
CORPORATION™

Choose Excellence, Expect Value

Surgical Specialties Corporation is the world's largest independent suture manufacturer. While headquartered in Westwood, Massachusetts, SSC is a global company with development and manufacturing facilities in the United States, Germany, Mexico, China, and the United Kingdom.

We design and manufacture high performance wound closure products including QUILL™, Sharpoint™, and LOOK™ brand sutures.

We are pleased to offer the new Sharpoint™ PLUS suture portfolio featuring German engineered needles.

Our philosophy is that performance doesn't have to be expensive.



Table of Contents

Suture Material Selection	3-4
Advanced Needle Technology	5
Needle Anatomy and Terminology	6
Needle Points and Cutting Edges	7
Needle Charts.....	8
 Absorbable Sutures	
Plain Gut.....	9-10
Chromic Gut	11-12
PolySyn FA™ (polyglycolic acid)	13
Monoderm™ (glycolide and ε-caprolactone)	14-15
PolySyn™ (polyglycolic acid).....	16-19
PDO (polydioxanone)	20-21
 Non-Absorbable Sutures	
Silk.....	22-23
Nylon™	24-25
Polypropylene.....	26-29
Polyviolene™ XL (polyester)- Coated.....	30-31
 Cross Reference Guide	
Ethicon Cross Reference.....	32-40
Covidien Cross Reference.....	41-45



SharpPoint ^{PLUS}
Surgical Sutures



Guide to Suture Materials

Absorbable

Material	Material Description	BSR Profile	Mass Absorption
GUT  Comparable to: Surgical Gut Suture - Plain	Twisted • Plain & Chromic • Packaged in Wetting Solution Common uses: General Closure, OB/Gyn, Urology, Bowel Anastomosis	In vivo strength retention: Varies	Varies depending upon type of suture (plain vs. chromic) and site of placement
CHROMIC GUT  Comparable to: Surgical Gut Suture - Chromic	Twisted • Plain & Chromic • Packaged in Wetting Solution Common uses: General Closure, OB/Gyn, Urology, Bowel Anastomosis	In vivo strength retention: Varies	Varies depending upon type of suture (plain vs. chromic) and site of placement
POLYSYN FA™ (polyglycolic acid)  • VICRYL RAPIDE™ * Velosorb™ Fast	Braided • Undyed (clear) • Coated Common uses: General Closure, OB/Gyn, Urology	In vivo strength retention: 48% at 7 days 0% at 14 days	Essentially complete between 49 and 63 days
MONODERM™ (polyglycolic acid and ε-caprolactone)  • MONOCRYL™ * Caprosyn™/Biosyn™	Monofilament • Undyed (clear) • Dyed (violet) Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 42% to 76% at 7 days 36% to 52% at 14 days	Essentially complete by 90 days
POLYSYN™ (polyglycolic acid)  • VICRYL™ * Polysorb™/ Dexon™	Braided • Undyed (clear) • Dyed (violet) Monofilament (9-0 & 10-0 only) • Undyed (clear) • Dyed (violet) • Coated & Uncoated Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 50% at 14 days 20% at 21 days	Essentially complete between 50 and 90 days
PDO (polydioxanone)  • PDS II™ * Maxon™	Monofilament • Dyed (violet) Common uses: Plastic Surgery, Orthopaedic, OB/Gyn, Gastrointestinal	In vivo strength retention: 80% to 90% at 14 days 60% to 82% at 28 days 47% to 79% at 42 days	Essentially complete between 180 and 220 days

• Comparable Ethicon™ Material * Comparable Covidien™ Material


Non-Absorbable

Material	Material Description	Qualities
SILK  • PERMA-HAND™ Silk * SofSilk™	Braided • Undyed (white) • Dyed (black) • Coated Common uses: General & Skin Closure, Gastrointestinal, Cardiovascular Surgery, Plastic Surgery, Ophthalmology	• Excellent handling and tying characteristics
NYLON™  • ETHILON™ * Monosof™/ Dermalon™	Monofilament • Undyed (clear) • Dyed (blue and black) • Packaged dry or in wetting solution Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery, Gastrointestinal	• Passes easily through tissue
POLYPROPYLENE  • PROLENE™ * Surgipro™/Surgipro™ II	Monofilament • Dyed (blue) Common uses: General & Skin Closure, Neurosurgery, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery	• Inert • Easy to handle • Passes easily through tissue
POLYVIOLENE™ (UNCOATED) (poly (ethylene-terephthalate))  • MERSILENE™ * Surgidac™	Braided • Undyed (white) • Dyed (green) • Uncoated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security
POLYVIOLENE™ XL (COATED)  • ETHIBOND™ * Tricon™	Braided • Undyed (white) • Dyed (green) • Coated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security

• Comparable Ethicon™ Material * Comparable Covidien™ Material


Advanced Needle Technology

Quill™ needle performance is an output of its design. Materials, geometry, and coatings work together to produce a high performance surgical instrument.




Proprietary Stainless Steel Alloy


Quill™ needles are forged from an alloy of Iron, Nickel, and Chromium. Each element serves a role in delivering strength, ductility, and hardness.



Strength



Ductility



Hardness

Needle Terminology

Length: Distance measured along the circumference of the needle, from the point to the swage. It is depicted on the packaging.

Shape: Different surgeries or surgical techniques demand different needle shapes.

Point: The sharpened part of the needle that first penetrates the tissue. The point consists of the area from the sharpened tip of the needle to the widest cross section of the needle body. The point determines ease of penetration and the initial size and shape of the hole.


Swage: Connection point of the suture and the needle. Micro teeth grip the suture material and prevent it from slipping out.

Wire Diameter: Gauge or thickness of the wire originally used to form the needle body.


Body: The part between the point and the swage of the needle.

Included Angle (Curvature): The fraction of the circle included in the curve of the needle. This is depicted on the packaging.


Needle Anatomy




Needle Shapes and General Uses




Straight (aka: Keith): Skin Closure



3/8 Circle: Skin Closure, Ophthalmic




1/2 Circle: Skin Closure, GI, Cardiovascular, Ophthalmic



5/8 Circle: Robotics and Laparoscopic

Needle Points And Cutting Edges

Needle Type	Point Geometry	Description	Commonly Used For
Precision Conventional Cutting		Triangular point with cutting edge on the inner curvature and precision-honed for added sharpness	General skin closure, subcutaneous tissue, sometimes for ophthalmic surgery, plastic or reconstructive surgery
Reverse Cutting		Triangular point with cutting edge on the outer curvature	Skin closure, ligamentous or fibrous tissues
Precision Reverse Cutting		Triangular point with cutting edge on the outer curvature and precision-honed for added sharpness	Primarily used in plastic surgery for delicate procedures and where a good cosmetic result is desired
Taper Point		Round shaft, straight or curved, taper point with no cutting edge	Soft tissue closure such as fascia, vascular, gastrointestinal, and most soft tissue below the skin
Taper Cutting		Four (4) cutting edges with a taper point	For penetration through tough tissue
Spatula		Spatula needle with the cutting edge on the inner or outer curvature	Ophthalmic and microsurgery



Precision Conventional Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
PC-1	PC-1	PC-13	
PC-3	PC-3	PC-11	
PC-5	PC-5	PC-12	

Precision Reverse Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
P-1	P-1	P-10	
G-1	G-1	G-1	
P-3	P-3	P-13	
G-3	G-3	G-3	
PS-3	PS-3	P-11	
PS-2	PS-2	P-12	
PS-1	PS-1	P-14	
PS	PS	P-15	
1/2 Circle			
PS-4	PS-4	P-24	

Reverse Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
FS-2	FS-2	C-13	
FS-1	FS-1	C-14	
FS	FS	C-15	
FSL	FSL	C-16	
FSLX	FSLX	C-17	
1/2 Circle			
CP-2	CP-2	GS-10	
OS-4	OS-4	HOS-10	
OS-6	OS-6	HOS-11	
CP-1	CP-1	GS-11	


Straight Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
SC-1	SC-1		
KS	KS	SC-2	

Taper Cutting Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
CC	CC	KV-1	
CC-1	CC-1	•	
V-4	V-4	KV-4	
1/2 Circle			
V-5	V-5	KV-5	
V-7	V-7	KV-7	
V-34	V-34	KV-34	

Taper Point Needles

Surgical Specialties	Ethicon™	Covidien™	
3/8 Circle			
BV-1	BV-1	CV-1	
BV	BV	CV	
C-1	C-1	CV-11	
BB	BB	CV-15	
1/2 Circle			
RB-2	RB-2	CV-22	
RB-1	RB-1	CV-23	
SH-1	SH-1	CV-25	
SH	SH	V-20	
CT-2	CT-2	GS-22	
MO-6	MO-6	HGS-22	
MH	MH	V-26	
CT-1	CT-1	GS-21	
CT	CT	GS-24X	
CTX	CTX	GS-25	
			
5/8 Circle			
UR-6	UR-6	GU-46	
UR-5	UR-5	GU-45	

Plain Gut

Precision Conventional Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-1  13mm 3/8 Circle	Express Gut	18"/45cm	B1916N	B1915N				

Precision Reverse Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
G-1  11mm 3/8 Circle	Plain Gut	18"/45cm	B774N					
G-1  11mm 3/8 Circle	Plain Gut	18"/45cm	B770N					
P-3  13mm 3/8 Circle	Plain Gut	18"/45cm		B686N	B1644N			
G-3  13mm 3/8 Circle	Plain Gut	18"/45cm		B772N				
PS-2  19mm 3/8 Circle	Plain Gut	18"/45cm			B1627N			

Reverse Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2  19mm 3/8 Circle	Plain Gut	27"/70cm			B821N	B822N		

N=box of 12; D=box of 36

Straight Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
SC-1  13mm Straight	Plain Gut	18"/45cm			B1824N			
SC-1  13mm Straight	Plain Gut	18"/45cm			B1828N			

Taper Point

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
CT-1  36mm 1/2 Circle	Plain Gut	27"/30cm				B842N	B843N	
CT  40mm 1/2 Circle	Plain Gut	27"/30cm				B852N	B853N	

Chromic Gut

Precision Reverse Cutting ▽				Size				
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
G-1  11mm 3/8 Circle	Chromic Gut	18"/45cm	C796N					
G-1  11mm 3/8 Circle	Chromic Gut	18"/45cm	C790N					
P-3  13mm 3/8 Circle	Chromic Gut	18"/45cm		C687N	C1654N			
G-3  13mm 3/8 Circle	Chromic Gut	18"/45cm		C792N	C793N			
PS-3  16mm 3/8 Circle	Chromic Gut	18"/45cm		C1636N				
PS-4  16mm 1/2 Circle	Chromic Gut	18"/45cm			C1643N			
PS-2  19mm 3/8 Circle	Chromic Gut	18"/45cm			C1637N			
		27"/70cm				C1638N		

Reverse Cutting ▽				Size				
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2  19mm 3/8 Circle	Chromic Gut	27"/70cm			C635N	C636N		

N=box of 12; D=box of 36

Taper Point 			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
RB-1  17mm 1/2 Circle	Chromic Gut	27"/70cm	C202N	C203N	C204N			
SH-1  22mm 1/2 Circle	Chromic Gut	27"/70cm		C181N	C182N			
SH  26mm 1/2 Circle	Chromic Gut	27"/70cm		C121N	C122N	C123N	C124N	
CT-1  36mm 1/2 Circle	Chromic Gut	27"/70cm			C810N	C811N	C812N	C813N
		36"/90cm			C922N	C923N	C924N	C925N
CT  40mm 1/2 Circle	Chromic Gut	27"/70cm					C802N	
		36"/90cm				C913N	C914N	C915N
CTX  48mm 1/2 Circle	Chromic Gut	36"/90cm					C904N	C905N

PolySyn FA™ / Polyglycolic Acid

Precision Reverse Cutting ▽

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-3  13mm 3/8 Circle	Undyed Braided	18"/45cm		M493N	M494N			
PS-2  19mm 3/8 Circle	Undyed Braided	18"/45cm			M496N	M497N		
		27"/70cm			M426N			

Taper Point ◎

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
CT-1  36mm 1/2 Circle	Undyed Braided	36"/90cm				M944N	M945N	

N=box of 12; D=box of 36

Precision Conventional Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-1  13mm 3/8 Circle	Undyed Monofilament	18"/45cm		Y834N				
PC-3  16mm 3/8 Circle	Undyed Monofilament	18"/45cm		Y844N	Y845N			
PC-5  19mm 3/8 Circle	Undyed Monofilament	18"/45cm			Y823N			

Precision Reverse Cutting

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-3  13mm 3/8 Circle	Undyed Monofilament	18"/45cm	Y492N	Y493N	Y494N			
	Violet Monofilament	18"/45cm		Y463N	Y464N			
PS-2  19mm 3/8 Circle	Undyed Monofilament	18"/45cm		Y495N	Y496N	Y497N		
		27"/70cm			Y426N	Y427N		
	Violet Monofilament	18"/45cm			Y513N			
PS-1  24mm 3/8 Circle	Undyed Monofilament	18"/45cm			Y682N			
		27"/70cm			Y935N	Y936N		

N=box of 12; D=box of 36

Monoderm™

Reverse Cutting ▽

			Size				
	Material	Length	5-0	4-0	3-0	2-0	0
CP-1  36mm 1/2 Circle	Undyed Monofilament	27"/70cm				Y266N	

Straight Cutting ▽

			Size			
	Material	Length	5-0	4-0	3-0	2-0
KS  60mm Straight	Undyed Monofilament	27"/70cm			Y523N	

Taper Point ●

			Size				
	Material	Length	4-0	3-0	2-0	0	1
RB-1  17mm 1/2 Circle	Undyed Monofilament	27"/70cm	Y214N	Y215N			
SH  26mm 1/2 Circle	Undyed Monofilament	27"/70cm	Y415N	Y416N	Y417N		
	Violet Monofilament	27"/70cm		Y316N			
CT-1  36mm 1/2 Circle	Undyed Monofilament	36"/90cm			Y945N	Y946N	
CTX  48mm 1/2 Circle	Violet Monofilament	36"/90cm				Y398N	Y399N

N=box of 12; D=box of 36

PolySyn™/Polyglycolic Acid

Precision Conventional Cutting

Precision Conventional Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-3 16mm 3/8 Circle 		Undyed Braided	18"/45cm			G845N			
PC-5 19mm 3/8 Circle 		Undyed Braided	18"/45cm				G824N		

Precision Reverse Cutting

Precision Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-1 11mm 3/8 Circle 		Undyed Braided	18"/45cm	G489N					
P-3 13mm 3/8 Circle 		Undyed Braided	18"/45cm	G492N	G493N	G494N			
		Violet Braided	18"/45cm		G463N	G464N			
PS-3 16mm 3/8 Circle 		Undyed Braided	18"/45cm		G500N				
PS-4 16mm 1/2 Circle 		Undyed Braided	18"/45cm			G507N			
PS-2 19mm 3/8 Circle 		Undyed Braided	18"/45cm		G495N	G496N	G497N		
			27"/70cm			G426N	G427N	G428N	
PS-1 24mm 3/8 Circle 		Undyed Braided	18"/45cm				G683N		
			27"/70cm			G935N	G936N		

N=box of 12; D=box of 36

PolySyn™/Polyglycolic Acid

Reverse Cutting ▽

			Size				
	Material	Length	4-0	3-0	2-0	0	1
FS-2  19mm 3/8 Circle	Undyed Braided	27"/70cm	G422N	G423N			
	Violet Braided	18"/45cm	G392N	G393N			
FS-1  24mm 3/8 Circle	Undyed Braided	27"/70cm	G441N	G442N	G443N		
CP-2  26mm 1/2 Circle	Undyed Braided	27"/70cm			G869N		
CP-1  36mm 1/2 Circle	Undyed Braided	27"/70cm			G266N	G267N	G268N
OS-6  36mm 1/2 Circle	Undyed Braided	27"/70cm					G535N

Straight Cutting ▽

			Size				
	Material	Length	6-0	5-0	4-0	3-0	2-0
KS  60mm Straight	Undyed Braided	27"/70cm			G662N		

Taper Cutting Ⓢ

			Size				
	Material	Length	5-0	4-0	3-0	2-0	0
V-34  36mm 1/2 Circle	Undyed Braided	36"/90cm				G917N	

N=box of 12; D=box of 36

Taper Point 				Size				
		Material	Length	4-0	3-0	2-0	0	1
RB-1 17mm 1/2 Circle 		Undyed Braided	27"/70cm	G214N	G215N			
SH 26mm 1/2 Circle 		Undyed Braided	27"/70cm	G415N	G416N	G417N		
		Violet Braided	27"/70cm		G316N	G317N		
CT-2 26mm 1/2 Circle 		Undyed Braided	27"/70cm			G269N	G270N	
		Violet Braided	27"/70cm			G333N	G334N	
UR-6 26mm 5/8 Circle 		Violet Braided	27"/70cm			G602N	G603N	
CT-1 36mm 1/2 Circle 		Undyed Braided	27"/70cm		G258N	G259N	G260N	G261N
			36"/90cm		G944N	G945N	G946N	G947N
		Violet Braided	27"/70cm		G338N	G339N	G340N	G341N
			36"/90cm		G344N	G345N	G346N	G347N
UR-5 36mm 5/8 Circle 		Violet Braided	27"/70cm				G376N	

PolySyn™/Polyglycolic Acid

Taper Point 			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
CT  40mm 1/2 Circle	Undyed Braided	36"/90cm				G957N	G958N	
	Violet Braided	36"/90cm					G358N	
CTX  48mm 1/2 Circle	Undyed Braided	36"/90cm					G978N	G977N
	Violet Braided	36"/90cm					G370N	G371N

N=box of 12; D=box of 36

Precision Reverse Cutting ▽

Precision Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-3 13mm 3/8 Circle 		Violet Monofilament	18"/45cm		R493N	R494N			
PS-2 19mm 3/8 Circle 		Violet Monofilament	18"/45cm			R513N/ R496N	R497N		

Reverse Cutting ▽

Reverse Cutting 				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
FS-2 19mm 3/8 Circle		Violet Monofilament	27"/70cm			R422N	R423N		
FS-1 24mm 3/8 Circle		Violet Monofilament	27"/70cm					R443N	

Taper Point ◎

Taper Point 			Size					
		Material	Length	7-0	6-0	5-0	4-0	3-0
BV-1 9.3mm 3/8 Circle 		Violet Monofilament	24"/60cm	R155N				
			30"/75cm	R135N	R117N			
C-1 13mm 3/8 Circle 		Violet Monofilament	30"/75cm		R127N	R126N		
RB-2 13mm 1/2 Circle 		Violet Monofilament	30"/75cm			R148N		
RB-1 17mm 1/2 Circle 		Violet Monofilament	27"/70cm			R303N	R304N	R305N

N=box of 12; D=box of 36

Taper Point

				Size					
		Material	Length	5-0	4-0	3-0	2-0	0	1
RB-1 17mm 1/2 Circle		Violet Monofilament	30"/75cm	R320N					
SH 26mm 1/2 Circle		Violet Monofilament	27"/70cm			R316N	R317N		
CT-2 26mm 1/2 Circle		Violet Monofilament	27"/70cm				R333N	R334N	
CT-1 36mm 1/2 Circle		Violet Monofilament	27"/70cm				R339N	R340N	R341N
			36"/90cm					R346N	R347N
CT 40mm 1/2 Circle		Violet Monofilament	36"/90cm					R358N	R359N
CTX 48mm 1/2 Circle		Violet Monofilament	36"/90cm					R370N	R371N
			60"/120cm					R990N	

N=box of 12; D=box of 36

Precision Reverse Cutting 			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
P-1 11mm 3/8 Circle 	Black Braided	18"/45cm	D639N					
G-1 11mm 3/8 Circle 	Black Braided	18"/45cm	D780N					
P-3 13mm 3/8 Circle 	Black Braided	18"/45cm		D640N	D641N			
G-3 13mm 3/8 Circle 	Black Braided	18"/45cm			D783N			
PS-1 24mm 3/8 Circle 	Black Braided	18"/45cm				D1684N		
PS 26mm 3/8 Circle 	Black Braided	18"/45cm					D1588N	
PSL 30mm 3/8 Circle 	Black Braided	18"/45cm						D580N
		30"/75cm						D590N

Reverse Cutting 			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
FS-1 24mm 3/8 Circle 	Black Braided	18"/45cm			D684N			
FS 26mm 3/8 Circle 	Black Braided	18"/45cm				D685N		
FSL 30mm 3/8 Circle 	Black Braided	18"/45cm					D678N	
		30"/75cm					D680N	

N=box of 12; D=box of 36

Silk

Straight Cutting ▽

			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
KS  60mm Straight	Black Braided	30"/75cm				D622N	D623N	

Taper Point ◎

			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
SH  26mm 1/2 Circle	Black Braided	30"/75cm			D832N	D833N	D834N	
CT-1  36mm 1/2 Circle	Black Braided	30"/75cm					D424N	

N=box of 12; D=box of 36

Precision Conventional Cutting

			Size				
Needle	Material	Length	6-0	5-0	4-0	3-0	2-0
PC-1  13mm 3/8 Circle	Black Monofilament	18"/45cm	A1856N	A1855N			
	Black Monofilament Pliable	18"/45cm		A1955N			
PC-3  16mm 3/8 Circle	Black Monofilament	18"/45cm		A1865N	A1864N		
	Black Monofilament Pliable	18"/45cm		A1965N			
PC-5  19mm 3/8 Circle	Black Monofilament	18"/45cm			A1894N	A1893N	
	Black Monofilament Pliable	18"/45cm			A1994N		

Precision Reverse Cutting

			Size				
Needle	Material	Length	6-0	5-0	4-0	3-0	2-0
P-1  11mm 3/8 Circle	Black Monofilament	18"/45cm	A697N				
P-3  13mm 3/8 Circle	Black Monofilament	18"/45cm	A1698N	A698N	A699N		
	Undyed Monofilament	18"/45cm		A690N	A691N		
PS-4  16mm 1/2 Circle	Black Monofilament	18"/45cm			A1662N		
PS-3  16mm 3/8 Circle	Black Monofilament	18"/45cm	A1665N	A1668N			
PS-2  19mm 3/8 Circle	Black Monofilament	18"/45cm		A1666N	A1667N	A1669N	A593N
	Blue Monofilament	18"/45cm			AG667N		

N=box of 12; D=box of 36

Nylon

Precision Reverse Cutting ▼

			Size				
	Material	Length	6-0	5-0	4-0	3-0	2-0
PS-1 24mm 3/8 Circle 	Black Monofilament	18"/45cm				A1663N	
PS 26mm 3/8 Circle 	Black Monofilament	18"/45cm					A585N
PSL 30mm 3/8 Circle 	Black Monofilament	30"/75cm				A1691N	
PSLX 36mm 3/8 Circle 	Black Monofilament	30"/75cm				A1683N	A1697N

Reverse Cutting ▼

			Size				
	Material	Length	5-0	4-0	3-0	2-0	0
FS-2 19mm 3/8 Circle 	Black Monofilament	18"/45cm	A661N	A662N			
FS-1 24mm 3/8 Circle 	Black Monofilament	18"/45cm		A1629N	A663N		
		30"/75cm			A669N		
FS 26mm 3/8 Circle 	Black Monofilament	18"/45cm				A664N	
FSL 30mm 3/8 Circle 	Black Monofilament	30"/75cm			A1671N		
FSLX 36mm 3/8 Circle 	Black Monofilament	30"/75cm			A1673N	A1674N	

N=box of 12; D=box of 36

Precision Conventional Cutting



Precision Conventional Cutting				Size					
		Material	Length	6-0	5-0	4-0	3-0	2-0	0
PC-1 13mm 3/8 Circle		Blue Monofilament	18"/45cm	J8617N	J8618N				
PC-3 16mm 3/8 Circle		Blue Monofilament	18"/45cm		J8635N	J8634N			
PC-5 19mm 3/8 Circle		Blue Monofilament	18"/45cm			J8631N	J8632N		

Precision Reverse Cutting



Precision Reverse Cutting				Size				
Material		Length	7-0	6-0	5-0	4-0	3-0	
P-1 11mm 3/8 Circle 	Blue Monofilament	18"/45cm		J8697N				
P-3 13mm 3/8 Circle 	Blue Monofilament	18"/45cm		J8695N	J8698N	J8699N		
PS-3 16mm 3/8 Circle 	Blue Monofilament	18"/45cm		J8680N	J8681N			
PS-2 19mm 3/8 Circle 	Blue Monofilament	18"/45cm			J8686N	J8682N	J8687N	
PS-1 24mm 3/8 Circle 	Blue Monofilament	18"/45cm					J8663N	

N=box of 12; D=box of 36

Polypropylene

Reverse Cutting ▽

Reverse Cutting 				Size				
		Material	Length	6-0	5-0	4-0	3-0	2-0
FS-2 19mm 3/8 Circle 		Blue Monofilament	18"/45cm		J8661N	J8683N	J8665N	
FS-1 24mm 3/8 Circle 		Blue Monofilament	18"/45cm				J8684N	
FS 26mm 3/8 Circle 		Blue Monofilament	18"/45cm					J8685N

Taper Point ⊙

Taper Point 			Size		
	Material	Length	7-0	6-0	5-0
BV-1  9.3mm 3/8 Circle	Blue Monofilament	18"/45cm	J8701D	J8806D	
		24"/60cm	J8702D/ J8704D	J8805D	
			J8747N	J8807D	
		30"/75cm	J8703D	J8709D	
			J8755N		
C-1  13mm 3/8 Circle	Blue Monofilament	18"/45cm		J8718D	J8717D
		24"/60cm		J8726D	J8725D
		30"/75cm		J8706D	
		36"/90cm			J8720D
BV  11mm 3/8 Circle	Blue Monofilament	24"/60cm		J8610N	
		30"/75cm		J8776N	


N=box of 12; D=box of 36

Taper Point 			Size					
	Material	Length	6-0	5-0	4-0	3-0	2-0	0
RB-2  13mm 1/2 Circle	Blue Monofilament	18"/45cm	J8714D	J8713D				
		24"/60cm		J8554D				
		30"/75cm	J8711D	J8710D				
BB  17mm 3/8 Circle	Blue Monofilament	36"/90cm		J8580D	J8581D			
RB-1  17mm 1/2 Circle	Blue Monofilament	18"/45cm		J8756D				
		24"/60cm		J8555D				
		36"/90cm		J8556D	J8557D	J8558D		
SH-1  22mm 1/2 Circle	Blue Monofilament	30"/75cm			J8526D	J8762D		
SH  26mm 1/2 Circle	Blue Monofilament	30"/75cm				J8832N	J8833N	J8834N
SH  26mm 1/2 Circle	Blue Monofilament	36"/90cm			J8521D	J8522D	J8523D	
		48"/122cm				J8534D	J8533D	
MO-6  26mm 1/2 Circle	Blue Monofilament	30"/75cm						J8418N
CT-2  26mm 1/2 Circle	Blue Monofilament	30"/75cm					J8411N	J8412N

N=box of 12; D=box of 36

Polypropylene

Taper Point

	Material	Length	3-0	2-0	0	1	2
MH  36mm 1/2 Circle	Blue Monofilament	36"/90cm	J8842D	J8843D			
		48"/122cm		J8853D			
CT-1  36mm 1/2 Circle	Blue Monofilament	30"/75cm		J8423N	J8424N	J8425N	
CT  40mm 1/2 Circle	Blue Monofilament	30"/75cm			J8434N	J8435N	
CTX  48mm 1/2 Circle	Blue Monofilament	30"/75cm				J8455N	

Taper Cutting

			Size					
	Material	Length	7-0	6-0	5-0	4-0	3-0	2-0
CC  9.3mm 3/8 Circle	Blue Monofilament	24"/60cm	J8704D	J8807D				
CC-1  13mm 3/8 Circle	Blue Monofilament	30"/75cm		J8707D				
V-5  17mm 1/2 Circle	Blue Monofilament	36"/90cm				J8935D		
V-7  26mm 1/2 Circle	Blue Monofilament	36"/90cm				J8975D	J8976D	J8977D

N=box of 12; D=box of 36

Polyviolene™ XL / Polyester - Coated

Reverse Cutting ▽

			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	2
OS-4  22mm 1/2 Circle	Green Braided	30"/75cm						E519N

Taper Point ◎


			Size					
	Material	Length	5-0	4-0	3-0	2-0	0	1
RB-1  17mm 1/2 Circle	Green Braided	30"/75cm		EX551D				
		36"/90cm	EX556D	EX557D	EX558D			
SH-1  22mm 1/2 Circle	Green Braided	36"/90cm			EX762D	EX763D		
SH  26mm 1/2 Circle	Green Braided	30"/75cm				EX563D		
		36"/90cm				EX523D	E524N	
	White Braided	36"/90cm				EX513D		
SH  26mm 1/2 Circle	Green Braided	30"/75cm				E833N	E834N	
CT-2  26mm 1/2 Circle	Green Braided	30"/75cm					E412N	
CT-1  36mm 1/2 Circle	Green Braided	30"/75cm					E424N	E425N

N=box of 12; D=box of 36

Polyviolene™ XL / Polyester - Coated

Taper Cutting



			Size					
Needle	Material	Length	6-0	5-0	4-0	3-0	2-0	0
V-4  17mm 3/8 Circle	Green Braided	30"/75cm					EX927D	
V-7  26mm 1/2 Circle	Green Braided	30"/75cm					EX977D	

N=box of 12; D=box of 36

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
1588H	D1588N	Silk 2-0 Black 18" PS Precision Reverse Cutting 3/8C 26mm
1627H	B1627N	Plain Gut 4-0 18" PS-2 Precision Reverse Cutting 3/8C 19mm
1629H	A1629N	Nylon 4-0 Black 18" FS-1 Reverse Cutting 3/8C 24mm
1636G	C1636N	Chromic Gut 5-0 18" PS-3 Precision Reverse Cutting 3/8C 16mm
1637G	C1637N	Chromic Gut 4-0 18" PS-2 Precision Reverse Cutting 3/8C 19mm
1638H	C1638N	Chromic Gut 3-0 27" PS-2 Precision Reverse Cutting 3/8C 19mm
1643G	C1643N	Chromic Gut 4-0 18" PS-4 Precision Reverse Cutting 1/2C 16mm
1644G	B1644N	Plain Gut 4-0 18" P-3 Precision Reverse Cutting 3/8C 13mm
1654G	C1654N	Chromic Gut 4-0 18" P-3 Precision Reverse Cutting 3/8C 13mm
1662G	A1662N	Nylon 4-0 Black 18" PS-4 Precision Reverse Cutting 1/2C 16mm
1663G	A1663N	Nylon 3-0 Black 18" PS-1 Precision Reverse Cutting 3/8C 24mm
1665G	A1665N	Nylon 6-0 Black 18" PS-3 Precision Reverse Cutting 3/8C 16mm
1666G	A1666N	Nylon 5-0 Black 18" PS-2 Precision Reverse Cutting 3/8C 19mm
1667G	A1667N	Nylon 4-0 Black 18" PS-2 Precision Reverse Cutting 3/8C 19mm
1668G	A1668N	Nylon 5-0 Black 18" PS-3 Precision Reverse Cutting 3/8C 16mm
1669H	A1669N	Nylon 3-0 Black 18" PS-2 Precision Reverse Cutting 3/8C 19mm
1671H	A1671N	Nylon 3-0 Black 30" FSL Reverse Cutting 3/8C 30mm
1673H	A1673N	Nylon 3-0 Black 30" FSLX Reverse Cutting 3/8C 36mm
1674H	A1674N	Nylon 2-0 Black 30" FSLX Reverse Cutting 3/8C 36mm
1683H	A1683N	Nylon 3-0 Black 30" PSLX Reverse Cutting 3/8C 36mm
1684G	D1684N	Silk 3-0 Black 18" PS-1 Precision Reverse Cutting 3/8C 24mm
1691H	A1691N	Nylon 3-0 Black 30" PSL Reverse Cutting 3/8C 30mm
1697H	A1697N	Nylon 2-0 Black 30" PSLX Reverse Cutting 3/8C 36mm
1698G	A1698N	Nylon 6-0 Black 18" P-3 Precision Reverse Cutting 3/8C 13mm
1824H	B1824N	Plain Gut 4-0 18" SC-1 Reverse Cutting Straight 13mm
1828H	B1828N	Plain Gut 4-0 18" SC-1 Reverse Cutting Straight 13mm Double Armed
1855G	A1855N	Nylon 5-0 Black 18" PC-1 Precision Conventional Cutting 3/8C 13mm
1856G	A1856N	Nylon 6-0 Black 18" PC-1 Precision Conventional Cutting 3/8C 13mm
1864G	A1864N	Nylon 4-0 Black 18" PC-3 Precision Conventional Cutting 3/8C 16mm
1865G	A1865N	Nylon 5-0 Black 18" PC-3 Precision Conventional Cutting 3/8C 16mm
1893G	A1893N	Nylon 3-0 Black 18" PC-5 Precision Conventional Cutting 3/8C 19mm
1894G	A1894N	Nylon 4-0 Black 18" PC-5 Precision Conventional Cutting 3/8C 19mm
1915G	B1915N	Express Gut 5-0 18" PC-1 Precision Conventional Cutting 3/8C 13mm
1916G	B1916N	Express Gut 6-0 18" PC-1 Precision Conventional Cutting 3/8C 13mm
1955G	A1955N	Nylon (pliable) 5-0 Black 18" PC-1 Precision Conventional Cutting 3/8C 13mm
1965G	A1965N	Nylon (pliable) 5-0 Black 18" PC-3 Precision Conventional Cutting 3/8C 16mm
1994G	A1994N	Nylon (Pliable) 4-0 Black 18" PC-5 Precision Conventional Cutting 3/8C 19mm

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
580H	D580N	Silk 0 Black 18" PSL Reverse Cutting 3/8C 30mm
585H	A585N	Nylon 2-0 Black 18" PS Precision Reverse Cutting 3/8C 26mm
590H	D590N	Silk 0 Black 30" PSL Reverse Cutting 3/8C 30mm
593H	A593N	Nylon 2-0 Black 18" PS-2 Precision Reverse Cutting 3/8C 19mm
623H	D623N	Silk 2-0 Black 30" GS60 Straight Straight Cutting 60mm
635H	C635N	Chromic Gut 4-0 27" FS-2 Reverse Cutting 3/8C 19mm
636H	C636N	Chromic Gut 3-0 27" FS-2 Reverse Cutting 3/8C 19mm
639G	D639N	Silk 6-0 Black 18" P-1 Precision Reverse Cutting 3/8C 11mm
640G	D640N	Silk 5-0 Black 18" P-3 Precision Reverse Cutting 3/8C 13mm
641G	D641N	Silk 4-0 Black 18" P-3 Precision Reverse Cutting 3/8C 13mm
661H	A661N	Nylon 5-0 Black 18" FS-2 Reverse Cutting 3/8C 19mm
662G	A662N	Nylon 4-0 Black 18" FS-2 Reverse Cutting 3/8C 19mm
663G	A663N	Nylon 3-0 Black 18" FS-1 Reverse Cutting 3/8C 24mm
664G	A664N	Nylon 2-0 Black 18" FS Reverse Cutting 3/8C 26mm
669H	A669N	Nylon 3-0 Black 30" FS-1 Reverse Cutting 3/8C 24mm
678G	D678N	Silk 0 Black 18" FSL Reverse Cutting 3/8C 30mm
680H	D680N	Silk 0 Black 30" FSL Reverse Cutting 3/8C 30mm
684G	D684N	Silk 3-0 Black 18" FS-1 Reverse Cutting 3/8C 24mm
685G	D685N	Silk 2-0 Black 18" FS Reverse Cutting 3/8C 26mm
686G	B686N	Plain Gut 5-0 18" P-3 Precision Reverse Cutting 3/8C 13mm
687G	C687N	Chromic Gut 5-0 18" P-3 Precision Reverse Cutting 3/8C 13mm
690G	A690N	Nylon 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
691G	A691N	Nylon 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
697H	A697N	Nylon 6-0 Black 18" P-1 Precision Reverse Cutting 3/8C 11mm
698H	A698N	Nylon 5-0 Black 18" P-3 Precision Reverse Cutting 3/8C 13mm
699G	A699N	Nylon 4-0 Black 18" P-3 Precision Reverse Cutting 3/8C 13mm
770G	B770N	Plain Gut 6-0 18" G-1 Precision Reverse Cutting 3/8C 11mm Double Armed
772G	B772N	Plain Gut 5-0 18" G-3 Precision Reverse Cutting 3/8C 13mm Double Armed
774G	B774N	Plain Gut 6-0 18" G-1 Precision Reverse Cutting 3/8C 11mm
780G	D780N	Silk 6-0 Black 18" G-1 Precision Reverse Cutting 3/8C 11mm Double Armed
783G	D783N	Silk 4-0 Black 18" G-3 Precision Reverse Cutting 3/8C 13mm Double Armed
790G	C790N	Chromic Gut 6-0 18" G-1 Precision Reverse Cutting 3/8C 11mm Double Armed
792G	C792N	Chromic Gut 5-0 18" G-3 Precision Reverse Cutting 3/8C 13mm Double Armed
793G	C793N	Chromic Gut 4-0 18" G-3 Precision Reverse Cutting 3/8C 13mm Double Armed
796G	C796N	Chromic Gut 6-0 18" G-1 Precision Reverse Cutting 3/8C 11mm
8521H	J8521D	Polypropylene 4-0 Blue 36" SH Taperpoint 1/2C 26mm Double Armed
8522H	J8522D	Polypropylene 3-0 Blue 36" SH Taperpoint 1/2C 26mm Double Armed
8523H	J8523D	Polypropylene 2-0 Blue 36" SH Taperpoint 1/2C 26mm Double Armed

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
8526H	J8526D	Polypropylene 4-0 Blue 30" SH-1 Taperpoint 1/2C 22mm Double Armed
8533H	J8533D	Polypropylene 2-0 Blue 48" SH Taperpoint 1/2C 26mm Double Armed
8534H	J8534D	Polypropylene 3-0 Blue 48" SH Taperpoint 1/2C 26mm Double Armed
8554H	J8554D	Polypropylene 5-0 Blue 24" RB-2 Taperpoint 1/2C 13mm Double Armed
8555H	J8555D	Polypropylene 5-0 Blue 24" RB-1 Taperpoint 1/2C 17mm Double Armed
8556H	J8556D	Polypropylene 5-0 Blue 36" RB-1 Taperpoint 1/2C 17mm Double Armed
8557H	J8557D	Polypropylene 4-0 Blue 36" RB-1 Taperpoint 1/2C 17mm Double Armed
8558H	J8558D	Polypropylene 3-0 Blue 36" RB-1 Taperpoint 1/2C 17mm Double Armed
8580H	J8580D	Polypropylene 5-0 Blue 36" BB Taperpoint 3/8C 17mm Double Armed
8581H	J8581D	Polypropylene 4-0 Blue 36" BB Taperpoint 3/8C 17mm Double Armed
8617G	J8617N	Polypropylene 6-0 Blue 18" PC-1 Precision Conventional Cutting 3/8C 13mm
8618G	J8618N	Polypropylene 5-0 Blue 18" PC-1 Precision Conventional Cutting 3/8C 13mm
8631G	J8631N	Polypropylene 4-0 Blue 18" PC-5 Precision Conventional Cutting 3/8C 19mm
8632G	J8632N	Polypropylene 3-0 Blue 18" PC-5 Precision Conventional Cutting 3/8C 19mm
8634G	J8634N	Polypropylene 4-0 Blue 18" PC-3 Precision Conventional Cutting 3/8C 16mm
8635G	J8635N	Polypropylene 5-0 Blue 18" PC-3 Precision Conventional Cutting 3/8C 16mm
8661G	J8661N	Polypropylene 5-0 Blue 18" FS-2 Reverse Cutting 3/8C 19mm
8663G	J8663N	Polypropylene 3-0 Blue 18" PS-1 Precision Reverse Cutting 3/8C 24mm
8665G	J8665N	Polypropylene 3-0 Blue 18" FS-2 Reverse Cutting 3/8C 19mm
8680G	J8680N	Polypropylene 6-0 Blue 18" PS-3 Precision Reverse Cutting 3/8C 16mm
8681G	J8681N	Polypropylene 5-0 Blue 18" PS-3 Precision Reverse Cutting 3/8C 16mm
8682G	J8682N	Polypropylene 4-0 Blue 18" PS-2 Precision Reverse Cutting 3/8C 19mm
8683G	J8683N	Polypropylene 4-0 Blue 18" FS-2 Reverse Cutting 3/8C 19mm
8684G	J8684N	Polypropylene 3-0 Blue 18" FS-1 Reverse Cutting 3/8C 24mm
8685H	J8685N	Polypropylene 2-0 Blue 18" FS Reverse Cutting 3/8C 26mm
8686G	J8686N	Polypropylene 5-0 Blue 18" PS-2 Precision Reverse Cutting 3/8C 19mm
8687H	J8687N	Polypropylene 3-0 Blue 18" PS-2 Precision Reverse Cutting 3/8C 19mm
8695G	J8695N	Polypropylene 6-0 Blue 18" P-3 Precision Reverse Cutting 3/8C 13mm
8697G	J8697N	Polypropylene 6-0 Blue 18" P-1 Precision Reverse Cutting 3/8C 11mm
8698G	J8698N	Polypropylene 5-0 Blue 18" P-3 Precision Reverse Cutting 3/8C 13mm
8699G	J8699N	Polypropylene 4-0 Blue 18" P-3 Precision Reverse Cutting 3/8C 13mm
8701H	J8701D	Polypropylene 7-0 Blue 18" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8702H	J8702D	Polypropylene 7-0 Blue 24" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8703H	J8703D	Polypropylene 7-0 Blue 30" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8704H	J8704D	Polypropylene 7-0 Blue 24" CC Tapercut 3/8C 9.3mm Double Armed
8706H	J8706D	Polypropylene 6-0 Blue 30" C-1 Taperpoint 3/8C 13mm Double Armed
8707H	J8707D	Polypropylene 6-0 Blue 30" CC-1 Tapercut Cardio 3/8C 13mm Double Armed
8709H	J8709D	Polypropylene 6-0 Blue 30" BV-1 Taperpoint 3/8C 9.3mm Double Armed

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
8710H	J8710D	Polypropylene 5-0 Blue 30" RB-2 Taperpoint 1/2C 13mm Double Armed
8711H	J8711D	Polypropylene 6-0 Blue 30" RB-2 Taperpoint 1/2C 13mm Double Armed
8713H	J8713D	Polypropylene 5-0 Blue 18" RB-2 Taperpoint 1/2C 13mm Double Armed
8714H	J8714D	Polypropylene 6-0 Blue 18" RB-2 Taperpoint 1/2C 13mm Double Armed
8717H	J8717D	Polypropylene 5-0 Blue 18" C-1 Taperpoint 3/8C 13mm Double Armed
8718H	J8718D	Polypropylene 6-0 Blue 18" C-1 Taperpoint 3/8C 13mm Double Armed
8720H	J8720D	Polypropylene 5-0 Blue 36" C-1 Taperpoint 3/8C 13mm Double Armed
8725H	J8725D	Polypropylene 5-0 Blue 24" C-1 Taperpoint 3/8C 13mm Double Armed
8726H	J8726D	Polypropylene 6-0 Blue 24" C-1 Taperpoint 3/8C 13mm Double Armed
8756H	J8756D	Polypropylene 5-0 Blue 18" RB-1 Taperpoint 1/2C 17mm Double Armed
8762H	J8762D	Polypropylene 3-0 Blue 30" SH-1 Taperpoint 1/2C 22mm Double Armed
8805H	J8805D	Polypropylene 6-0 Blue 24" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8806H	J8806D	Polypropylene 6-0 Blue 18" BV-1 Taperpoint 3/8C 9.3mm Double Armed
8807H	J8807D	Polypropylene 6-0 Blue 24" CC Tapercut 3/8C 9.3mm Double Armed
8842H	J8842D	Polypropylene 3-0 Blue 36" MH Taperpoint 1/2C 36mm Double Armed
8843H	J8843D	Polypropylene 2-0 Blue 36" MH Taperpoint 1/2C 36mm Double Armed
8853H	J8853D	Polypropylene 2-0 Blue 48" MH Taperpoint 1/2C 36mm Double Armed
8935H	J8935D	Polypropylene 4-0 Blue 36" V-5 Tapercut Cardio 1/2C 17mm Double Armed
8975H	J8975D	Polypropylene 4-0 Blue 36" V-7 Tapercut Cardio 1/2C 26mm Double Armed
8976H	J8976D	Polypropylene 3-0 Blue 36" V-7 Tapercut Cardio 1/2C 26mm Double Armed
8977H	J8977D	Polypropylene 2-0 Blue 36" V-7 Tapercut Cardio 1/2C 26mm Double Armed
G182H	C182N	Chromic Gut 3-0 30" HR22 1/2 Circle Taper Point 22mm
G667G	AG667N	Nylon 4-0 Blue 18" PS-2 Precision Reverse Cutting 3/8C 19mm
H821H	B821N	Plain Gut 4-0 27" FS-2 Reverse Cutting 3/8C 19mm
H822H	B822N	Plain Gut 3-0 27" FS-2 Reverse Cutting 3/8C 19mm
J232H	G416N	PolySyn™ 3-0 Undyed 27" HR26 1/2 Circle Taper Point 26mm
J260H	G260N	PolySyn™ 0 Undyed 27" HR36 1/2 Circle Taper Point 36mm
J266H	G266N	PolySyn 2-0 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
J267H	G267N	PolySyn 0 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
J268H	G268N	PolySyn 1 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
J269H	G269N	PolySyn™ 2-0 Undyed 27" HR26 1/2 Circle Taper Point 26mm
J316H	G316N	PolySyn™ 3-0 Violet 27" HR26 1/2 Circle Taper Point 26mm
J317H	G317N	PolySyn™ 2-0 Violet 27" HR26 1/2 Circle Taper Point 26mm
J333H	G333N	PolySyn™ 2-0 Violet 27" HR26 1/2 Circle Taper Point 26mm
J339H	G339N	PolySyn™ 2-0 Violet 27" HR36 1/2 Circle Taper Point 36mm
J340H	G340N	PolySyn™ 0 Violet 27" HRX36 1/2 Circle Taper Point 36.6mm
J341H	G341N	PolySyn™ 1 Violet 27" HRX36 1/2 Circle Taper Point 36.6mm
J358H	G358N	PolySyn™ 0 Violet 36" HRX40 1/2 Circle Taper Point 40mm

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
J392H	G392N	PolySyn 4-0 Violet 18" FS-2 Reverse Cutting 3/8C 19mm
J393H	G393N	PolySyn 3-0 Violet 18" FS-2 Reverse Cutting 3/8C 19mm
J422H	G422N	PolySyn 4-0 Undyed 27" FS-2 Reverse Cutting 3/8C 19mm
J423H	G423N	PolySyn 3-0 Undyed 27" FS-2 Reverse Cutting 3/8C 19mm
J426H	G426N	PolySyn 4-0 Undyed 27" PS-2 Precision Reverse Cutting 3/8C 19mm
J427H	G427N	PolySyn 3-0 Undyed 27" PS-2 Precision Reverse Cutting 3/8C 19mm
J428H	G428N	PolySyn 2-0 Undyed 27" PS-2 Precision Reverse Cutting 3/8C 19mm
J441H	G441N	PolySyn 4-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
J442H	G442N	PolySyn 3-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
J443H	G443N	PolySyn 2-0 Undyed 27" FS-1 Reverse Cutting 3/8C 24mm
J463G	G463N	PolySyn 5-0 Violet 18" P-3 Precision Reverse Cutting 3/8C 13mm
J464G	G464N	PolySyn 4-0 Violet 18" P-3 Precision Reverse Cutting 3/8C 13mm
J489G	G489N	PolySyn 6-0 Undyed 18" P-1 Precision Reverse Cutting 3/8C 11mm
J492G	G492N	PolySyn 6-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
J493G	G493N	PolySyn 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
J494G	G494N	PolySyn 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
J495G	G495N	PolySyn 5-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
J496G	G496N	PolySyn 4-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
J497G	G497N	PolySyn 3-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
J500G	G500N	PolySyn 5-0 Undyed 18" PS-3 Precision Reverse Cutting 3/8C 16mm
J507G	G507N	PolySyn 4-0 Undyed 18" PS-4 Precision Reverse Cutting 1/2C 16mm
J535H	G535N	PolySyn 1 Undyed 27" OS-6 Reverse Cutting 1/2C 36mm
J683G	G683N	PolySyn 3-0 Undyed 18" PS-1 Precision Reverse Cutting 3/8C 24mm
J824G	G824N	PolySyn 3-0 Undyed 18" PC-5 Precision Conventional Cutting 3/8C 19mm
J845G	G845N	PolySyn 4-0 Undyed 18" PC-3 Precision Conventional Cutting 3/8C 16mm
J869H	G869N	PolySyn 2-0 Undyed 27" CP-2 Reverse Cutting 1/2C 26mm
J935H	G935N	PolySyn 4-0 Undyed 27" PS-1 Precision Reverse Cutting 3/8C 24mm
J936H	G936N	PolySyn 3-0 Undyed 27" PS-1 Precision Reverse Cutting 3/8C 24mm
J945H	G945N	PolySyn™ 2-0 Undyed 36" HR36 1/2 Circle Taper Point 36mm
K832H	D832N	Silk 3-0 Black 30" HR26 1/2 Circle Taper Point 26mm
VR426	M426N	PolySyn FA 4-0 Undyed 27" PS-2 Precision Reverse Cutting 3/8C 19mm
VR493	M493N	PolySyn FA 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
VR494	M494N	PolySyn FA 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
VR496	M496N	PolySyn FA 4-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
VR497	M497N	PolySyn FA 3-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
X519H	E519N	Polyviolene 2 Green 30" OS-4 Reverse Cutting 1/2C 22mm
Y266H	Y266N	Monoderm 2-0 Undyed 27" CP-1 Reverse Cutting 1/2C 36mm
Y426H	Y426N	Monoderm 4-0 Undyed 27" PS-2 Precision Reverse Cutting 3/8C 19mm

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
Y427H	Y427N	Monoderm 3-0 Undyed 27" PS-2 Precision Reverse Cutting 3/8C 19mm
Y463G	Y463N	Monoderm 5-0 Violet 18" P-3 Precision Reverse Cutting 3/8C 13mm
Y464G	Y464N	Monoderm 4-0 Violet 18" P-3 Precision Reverse Cutting 3/8C 13mm
Y492G	Y492N	Monoderm 6-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
Y493G	Y493N	Monoderm 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
Y494G	Y494N	Monoderm 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8C 13mm
Y495G	Y495N	Monoderm 5-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
Y496G	Y496N	Monoderm 4-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
Y497G	Y497N	Monoderm 3-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8C 19mm
Y513G	Y513N	Monoderm 4-0 Violet 18" PS-2 Precision Reverse Cutting 3/8C 19mm
Y682H	Y682N	Monoderm 4-0 Undyed 18" PS-1 Precision Reverse Cutting 3/8C 24mm
Y823G	Y823N	Monoderm 4-0 Undyed 18" PC-5 Precision Conventional Cutting 3/8C 19mm
Y834G	Y834N	Monoderm 5-0 Undyed 18" PC-1 Precision Conventional Cutting 3/8C 13mm
Y844G	Y844N	Monoderm 5-0 Undyed 18" PC-3 Precision Conventional Cutting 3/8C 16mm
Y845G	Y845N	Monoderm 4-0 Undyed 18" PC-3 Precision Conventional Cutting 3/8C 16mm
Y935H	Y935N	Monoderm 4-0 Undyed 27" PS-1 Precision Reverse Cutting 3/8C 24mm
Y936H	Y936N	Monoderm 3-0 Undyed 27" PS-1 Precision Reverse Cutting 3/8C 24mm
Z422H	R422N	PDO 4-0 Violet 27" FS-2 Reverse Cutting 3/8C 19mm
Z423H	R423N	PDO 3-0 Violet 27" FS-2 Reverse Cutting 3/8C 19mm
Z443H	R443N	PDO 2-0 Violet 27" FS-1 Reverse Cutting 3/8C 24mm
Z493G	R493N	PDO 5-0 Violet 18" P-3 Precision Reverse Cutting 3/8C 13mm
Z494G	R494N	PDO 4-0 Violet 18" P-3 Precision Reverse Cutting 3/8C 13mm
Z496G	R496N	PDO 4-0 Violet 18" PS-2 Precision Reverse Cutting 3/8C 19mm
Z497G	R497N	PDO 3-0 Violet 18" PS-2 Precision Reverse Cutting 3/8C 19mm
Z513G	R513N	PDO 4-0 Violet 18" PS-2 Precision Reverse Cutting 3/8C 19mm
842H	B842N	Plain Gut 3-0 27" CT-1 Taperpoint 1/2 Circle 36mm
843H	B843N	Plain Gut 2-0 27" CT-1 Taperpoint 1/2 Circle 36mm
852H	B852N	Plain Gut 3-0 27" CT Taperpoint 1/2 Circle 40mm
853H	B853N	Plain Gut 2-0 27" CT Taperpoint 1/2 Circle 40mm
G121H	C121N	Chromic Gut 4-0 27" SH Taperpoint 1/2 Circle 26mm
G122H	C122N	Chromic Gut 3-0 27" SH Taperpoint 1/2 Circle 26mm
G123H	C123N	Chromic Gut 2-0 27" SH Taperpoint 1/2 Circle 26mm
G124H	C124N	Chromic Gut 0 27" SH Taperpoint 1/2 Circle 26mm
G181H	C181N	Chromic Gut 4-0 27" SH-1 Taperpoint 1/2 Circle 22mm
U202H	C202N	Chromic Gut 5-0 27" RB-1 Taperpoint 1/2 Circle 17mm
U203H	C203N	Chromic Gut 4-0 27" RB-1 Taperpoint 1/2 Circle 17mm
U204H	C204N	Chromic Gut 3-0 27" RB-1 Taperpoint 1/2 Circle 17mm
802H	C802N	Chromic Gut 0 27" CT Taperpoint 1/2 Circle 40mm

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
810H	C810N	Chromic Gut 3-0 27" CT-1 Taperpoint 1/2 Circle 36mm
811H	C811N	Chromic Gut 2-0 27" CT-1 Taperpoint 1/2 Circle 36mm
812H	C812N	Chromic Gut 0 27" CT-1 Taperpoint 1/2 Circle 36mm
813H	C813N	Chromic Gut 1 27" CT-1 Taperpoint 1/2 Circle 36mm
904H	C904N	Chromic Gut 0 36" CTX Taperpoint 1/2 Circle 48mm
905H	C905N	Chromic Gut 1 36" CTX Taperpoint 1/2 Circle 48mm
913H	C913N	Chromic Gut 2-0 36" CT Taperpoint 1/2 Circle 40mm
914H	C914N	Chromic Gut 0 36" CT Taperpoint 1/2 Circle 40mm
915H	C915N	Chromic Gut 1 36" CT Taperpoint 1/2 Circle 40mm
922H	C922N	Chromic Gut 3-0 36" CT-1 Taperpoint 1/2 Circle 36mm
923H	C923N	Chromic Gut 2-0 36" CT-1 Taperpoint 1/2 Circle 36mm
924H	C924N	Chromic Gut 0 36" CT-1 Taperpoint 1/2 Circle 36mm
925H	C925N	Chromic Gut 1 36" CT-1 Taperpoint 1/2 Circle 36mm
424H	D424N	Silk 0 Black 30" CT-1 Taperpoint 1/2 Circle 36mm
622H	D622N	Silk 3-0 Black 30" KS Reverse Cutting Straight 60mm
K833H	D833N	Silk 2-0 Black 30" SH Taperpoint 1/2 Circle 26mm
K834H	D834N	Silk 0 Black 30" SH Taperpoint 1/2 Circle 26mm
X834H	E412N	Polyviolene 0 Green 30" CT-2 Taperpoint 1/2 Circle 26mm
X424H	E424N	Polyviolene 0 Green 30" CT-1 Taperpoint 1/2 Circle 36mm
X425H	E425N	Polyviolene 1 Green 30" CT-1 Taperpoint 1/2 Circle 36mm
X524H	E524N	Polyviolene 0 Green 36" SH Taperpoint 1/2 Circle 26mm Double Armed
X833H	E833N	Polyviolene 2-0 Green 30" SH Taperpoint 1/2 Circle 26mm
X834H	E834N	Polyviolene 0 Green 30" SH Taperpoint 1/2 Circle 26mm
J214H	G214N	PolySyn 4-0 Undyed 27" RB-1 Taperpoint 1/2 Circle 17mm
J215H	G215N	PolySyn 3-0 Undyed 27" RB-1 Taperpoint 1/2 Circle 17mm
J258H	G258N	PolySyn 3-0 Undyed 27" CT-1 Taperpoint 1/2 Circle 36mm
J259H	G259N	PolySyn 2-0 Undyed 27" CT-1 Taperpoint 1/2 Circle 36mm
J261H	G261N	PolySyn 1 Undyed 27" CT-1 Taperpoint 1/2 Circle 36mm
J270H	G270N	PolySyn 0 Undyed 27" CT-2 Taperpoint 1/2 Circle 26mm
J334H	G334N	PolySyn 0 Violet 27" CT-2 Taperpoint 1/2 Circle 26mm
J338H	G338N	PolySyn 3-0 Violet 27" CT-1 Taperpoint 1/2 Circle 36mm
J344H	G344N	PolySyn 3-0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
J345H	G345N	PolySyn 2-0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
J346H	G346N	PolySyn 0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
J347H	G347N	PolySyn 1 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
J370H	G370N	PolySyn 0 Violet 36" CTX Taperpoint 1/2 Circle 48mm
J371H	G371N	PolySyn 1 Violet 36" CTX Taperpoint 1/2 Circle 48mm
J376H	G376N	PolySyn 0 Violet 27" UR-5 Taperpoint 5/8C 36mm

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
J415H	G415N	PolySyn 4-0 Undyed 27" SH Taperpoint 1/2 Circle 26mm
J269H	G417N	PolySyn 2-0 Undyed 27" SH Taperpoint 1/2 Circle 26mm
J602H	G602N	PolySyn 2-0 Violet 27" UR-6 Taperpoint 5/8C 26mm
J603H	G603N	PolySyn 0 Violet 27" UR-6 Taperpoint 5/8C 26mm
J662H	G662N	PolySyn 4-0 Undyed 27" KS Reverse Cutting Straight 60mm
J917H	G917N	PolySyn 2-0 Undyed 36" V-34 Tapercut 1/2 Circle 36mm
J944H	G944N	PolySyn 3-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
J946H	G946N	PolySyn 0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
J947H	G947N	PolySyn 1 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
J957H	G957N	PolySyn 2-0 Undyed 36" CT Taperpoint 1/2 Circle 40mm
J958H	G958N	PolySyn 0 Undyed 36" CT Taperpoint 1/2 Circle 40mm
J977H	G977N	PolySyn 1 Undyed 36" CTX Taperpoint 1/2 Circle 48mm
J978H	G978N	PolySyn 0 Undyed 36" CTX Taperpoint 1/2 Circle 48mm
8833H	J8411N	Polypropylene 2-0 Blue 30" CT-2 Taperpoint 1/2 Circle 26mm
8834H	J8412N	Polypropylene 0 Blue 30" CT-2 Taperpoint 1/2 Circle 26mm
8834H	J8418N	Polypropylene 0 Blue 30" MO-6 Taperpoint 1/2 Circle 26mm
8423H	J8423N	Polypropylene 2-0 Blue 30" CT-1 Taperpoint 1/2 Circle 36mm
8424H	J8424N	Polypropylene 0 Blue 30" CT-1 Taperpoint 1/2 Circle 36mm
8425H	J8425N	Polypropylene 1 Blue 30" CT-1 Taperpoint 1/2 Circle 36mm
8434H	J8434N	Polypropylene 0 Blue 30" CT Taperpoint 1/2 Circle 40mm
8435H	J8435N	Polypropylene 1 Blue 30" CT Taperpoint 1/2 Circle 40mm
8455H	J8455N	Polypropylene 1 Blue 30" CTX Taperpoint 1/2 Circle 48mm
EP8610H	J8610N	Polypropylene 6-0 Blue 24" BV Taperpoint 3/8 Circle 11mm Double Armed
8702H	J8747N	Polypropylene 7-0 Blue 24" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
8755H	J8755N	Polypropylene 7-0 Blue 30" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
EP8776H	J8776N	Polypropylene 6-0 Blue 30" BV Taperpoint 3/8 Circle 11mm Double Armed
8410H	J8832N	Polypropylene 3-0 Blue 30" SH Taperpoint 1/2 Circle 26mm
8411H	J8833N	Polypropylene 2-0 Blue 30" SH Taperpoint 1/2 Circle 26mm
8834H	J8834N	Polypropylene 0 Blue 30" SH Taperpoint 1/2 Circle 26mm
VCP944H	M944N	PolySyn FA 3-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
VR945	M945N	PolySyn FA 2-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
Z117H	R117N	PDO 6-0 Violet 30" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
Z126H	R126N	PDO 5-0 Violet 30" C-1 Taperpoint 3/8 Circle 13mm Double Armed
Z127H	R127N	PDO 6-0 Violet 30" C-1 Taperpoint 3/8 Circle 13mm Double Armed
Z135H	R135N	PDO 7-0 Violet 30" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
Z148H	R148N	PDO 5-0 Violet 30" RB-2 Taperpoint 1/2 Circle 13mm Double Armed
Z155H	R155N	PDO 7-0 Violet 24" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
Z259H	R259N	PDO 2-0 Undyed 27" CT-1 Taperpoint 1/2 Circle 36mm

Ethicon Cross Reference Guide

Ethicon	SSC	Product Description
Z303H	R303N	PDO 5-0 Violet 27" RB-1 Taperpoint 1/2 Circle 17mm
Z304H	R304N	PDO 4-0 Violet 27" RB-1 Taperpoint 1/2 Circle 17mm
Z305H	R305N	PDO 3-0 Violet 27" RB-1 Taperpoint 1/2 Circle 17mm
Z316H	R316N	PDO 3-0 Violet 27" SH Taperpoint 1/2 Circle 26mm
Z333H	R317N	PDO 2-0 Violet 27" SH Taperpoint 1/2 Circle 26mm
Z320H	R320N	PDO 5-0 Violet 30" RB-1 Taperpoint 1/2 Circle 17mm Double Armed
Z333H	R333N	PDO 2-0 Violet 27" CT-2 Taperpoint 1/2 Circle 26mm
Z334H	R334N	PDO 0 Violet 27" CT-2 Taperpoint 1/2 Circle 26mm
Z339H	R339N	PDO 2-0 Violet 27" CT-1 Taperpoint 1/2 Circle 36mm
Z340H	R340N	PDO 0 Violet 27" CT-1 Taperpoint 1/2 Circle 36mm
Z341H	R341N	PDO 1 Violet 27" CT-1 Taperpoint 1/2 Circle 36mm
Z346H	R346N	PDO 0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
Z347H	R347N	PDO 1 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
Z358T	R358N	PDO 0 Violet 36" CT Taperpoint 1/2 Circle 40mm
Z359T	R359N	PDO 1 Violet 36" CT Taperpoint 1/2 Circle 40mm
Z370T	R370N	PDO 0 Violet 36" CTX Taperpoint 1/2 Circle 48mm
Z371T	R371N	PDO 1 Violet 36" CTX Taperpoint 1/2 Circle 48mm
Z416H	R416N	PDO 3-0 Undyed 27" SH Taperpoint 1/2 Circle 26mm
Z990G	R990N	PDO 0 Violet 60" CTX Taperpoint 1/2 Circle 48mm
Y214H	Y214N	Monoderm 4-0 Undyed 27" RB-1 Taperpoint 1/2 Circle 17mm
Y215H	Y215N	Monoderm 3-0 Undyed 27" RB-1 Taperpoint 1/2 Circle 17mm
Y316H	Y316N	Monoderm 3-0 Violet 27" SH Taperpoint 1/2 Circle 26mm
Y398H	Y398N	Monoderm 0 Violet 36" CTX Taperpoint 1/2 Circle 48mm
Y399H	Y399N	Monoderm 1 Violet 36" CTX Taperpoint 1/2 Circle 48mm
Y415H	Y415N	Monoderm 4-0 Undyed 27" SH Taperpoint 1/2 Circle 26mm
Y416H	Y416N	Monoderm 3-0 Undyed 27" SH Taperpoint 1/2 Circle 26mm
Y417H	Y417N	Monoderm 2-0 Undyed 27" SH Taperpoint 1/2 Circle 26mm
Y523H	Y523N	Monoderm 3-0 Undyed 27" KS Reverse Cutting Straight 60mm
Y945H	Y945N	Monoderm 2-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
Y946H	Y946N	Monoderm 0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm

Covidien Cross Reference Guide

Covidien	SSC	Product Description
8886291451	EX513D	Polyviolene XL 2-0 36" SH SH Taper Point 1/2 Circle 26mm
8886295351	EX513D	Polyviolene XL 2-0 36" SH SH Taper Point 1/2 Circle 26mm
CC923	Y945N	Monoderm 2-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
CC924	Y946N	Monoderm 0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
CC975	Y399N	Monoderm 1 Violet 36" CTX Taperpoint 1/2 Circle 48mm
CL904	G370N	PolySyn 0 Violet 36" CTX Taperpoint 1/2 Circle 48mm
CL905	G371N	PolySyn 1 Violet 36" CTX Taperpoint 1/2 Circle 48mm
CL914	G358N	PolySyn 0 Violet 36" CT Taperpoint 1/2 Circle 40mm
CL922	G344N	PolySyn 3-0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
CL923	G345N	PolySyn 2-0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
CL924	G346N	PolySyn 0 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
CL925	G347N	PolySyn 1 Violet 36" CT-1 Taperpoint 1/2 Circle 36mm
CL952	G944N	PolySyn 3-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
CL953	G945N	PolySyn 2-0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
CL954	G946N	PolySyn 0 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
CL955	G947N	PolySyn 1 Undyed 36" CT-1 Taperpoint 1/2 Circle 36mm
CL963	G957N	PolySyn 2-0 Undyed 36" CT Taperpoint 1/2 Circle 40mm
CL964	G958N	PolySyn 0 Undyed 36" CT Taperpoint 1/2 Circle 40mm
CL974	G978N	PolySyn 0 Undyed 36" CTX Taperpoint 1/2 Circle 48mm
CL975	G977N	PolySyn 1 Undyed 36" CTX Taperpoint 1/2 Circle 48mm
CM904	Y398N	Monoderm 0 Violet 36" CTX Taperpoint 1/2 Circle 48mm
CP411	J8411N	Polypropylene 2-0 Blue 30" CT-2 Taperpoint 1/2 Circle 26mm
CP412	J8412N	Polypropylene 0 Blue 30" CT-2 Taperpoint 1/2 Circle 26mm
CP412	J8418N	Polypropylene 0 Blue 30" MO-6 Taperpoint 1/2 Circle 26mm
CP423	J8423N	Polypropylene 2-0 Blue 30" CT-1 Taperpoint 1/2 Circle 36mm
CP424	J8424N	Polypropylene 0 Blue 30" CT-1 Taperpoint 1/2 Circle 36mm
CP425	J8425N	Polypropylene 1 Blue 30" CT-1 Taperpoint 1/2 Circle 36mm
CP434	J8434N	Polypropylene 0 Blue 30" CT Taperpoint 1/2 Circle 40mm
CP435	J8435N	Polypropylene 1 Blue 30" CT Taperpoint 1/2 Circle 40mm
CP455	J8455N	Polypropylene 1 Blue 30" CTX Taperpoint 1/2 Circle 48mm
CS424	D424N	Silk 0 Black 30" CT-1 Taperpoint 1/2 Circle 36mm
GS65M	D834N	Silk 0 Black 30" SH Taperpoint 1/2 Circle 26mm
GS66M	D833N	Silk 2-0 Black 30" SH Taperpoint 1/2 Circle 26mm
GS67M	D832N	Silk 3-0 Black 30" SH Taperpoint 1/2 Circle 26mm
SC5587G	Y463N	Monoderm 5-0 Violet 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SC5590G	Y464N	Monoderm 4-0 Violet 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
SC5616G	Y495N	Monoderm 5-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SC5617G	Y496N	Monoderm 4-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SC5618G	Y497N	Monoderm 3-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SC5627G	Y513N	Monoderm 4-0 Violet 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SC5688G	Y492N	Monoderm 6-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SC5689G	Y493N	Monoderm 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SC5690G	Y494N	Monoderm 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SDN5666G	AG667N	Nylon 4-0 Blue 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SL1625G	G489N	PolySyn 6-0 Undyed 18" P-1 Precision Reverse Cutting 3/8 Circle 11mm
SL1654	G824N	PolySyn 3-0 Undyed 18" PC-5 Precision Conventional Cutting 3/8 Circle 19mm
SL1657	G845N	PolySyn 4-0 Undyed 18" PC-3 Precision Conventional Cutting 3/8 Circle 16mm
SL5626	G500N	PolySyn 5-0 Undyed 18" PS-3 Precision Reverse Cutting 3/8 Circle 16mm
SL5627G	G496N	PolySyn 4-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SL5628G	G497N	PolySyn 3-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SL5633G	G507N	PolySyn 4-0 Undyed 18" PS-4 Precision Reverse Cutting 1/2 Circle 16mm
SL5636G	G495N	PolySyn 5-0 Undyed 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SL5640G	G683N	PolySyn 3-0 Undyed 18" PS-1 Precision Reverse Cutting 3/8 Circle 24mm
SL5687G	G493N	PolySyn 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SL5688G	G492N	PolySyn 6-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SL5691G	G494N	PolySyn 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SN1693	A1893N	Nylon 3-0 Black 18" PC-5 Precision Conventional Cutting 3/8 Circle 19mm
SN1964	A1864N	Nylon 4-0 Black 18" PC-3 Precision Conventional Cutting 3/8 Circle 16mm
SN1994	A1894N	Nylon 4-0 Black 18" PC-5 Precision Conventional Cutting 3/8 Circle 19mm
SN1994	A1994N	Nylon (Pliable) 4-0 Black 18" PC-5 Precision Conventional Cutting 3/8 Circle 19mm
SN3697	A697N	Nylon 6-0 Black 18" P-1 Precision Reverse Cutting 3/8 Circle 11mm
SN3965	A1865N	Nylon 5-0 Black 18" PC-3 Precision Conventional Cutting 3/8 Circle 16mm
SN3965	A1965N	Nylon (Pliable) 5-0 Black 18" PC-3 Precision Conventional Cutting 3/8 Circle 16mm
SN5662G	A1662N	Nylon 4-0 Black 18" PS-4 Precision Reverse Cutting 1/2 Circle 16mm
SN5663G	A1663N	Nylon 3-0 Black 18" PS-1 Precision Reverse Cutting 3/8 Circle 24mm
SN5666G	A1666N	Nylon 5-0 Black 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SN5667G	A1667N	Nylon 4-0 Black 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SN5668G	A1668N	Nylon 5-0 Black 18" PS-3 Precision Reverse Cutting 3/8 Circle 16mm
SN5669G	A1669N	Nylon 3-0 Black 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SN5670	A593N	Nylon 2-0 Black 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SN5690G	A690N	Nylon 5-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SN5691	A691N	Nylon 4-0 Undyed 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
SN5696G	A1698N	Nylon 6-0 Black 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SN5698G	A698N	Nylon 5-0 Black 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SN5699G	A699N	Nylon 4-0 Black 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SN658	A1629N	Nylon 4-0 Black 18" FS-1 Reverse Cutting 3/8 Circle 24mm
SN661G	A661N	Nylon 5-0 Black 18" FS-2 Reverse Cutting 3/8 Circle 19mm
SN662G	A662N	Nylon 4-0 Black 18" FS-2 Reverse Cutting 3/8 Circle 19mm
SN663G	A663N	Nylon 3-0 Black 18" FS-1 Reverse Cutting 3/8 Circle 24mm
SN669	A669N	Nylon 3-0 Black 30" FS-1 Reverse Cutting 3/8 Circle 24mm
SN671G	A1671N	Nylon 3-0 Black 30" FSL Reverse Cutting 3/8 Circle 30mm
SN671G	A1691N	Nylon 3-0 Black 30" PSL Reverse Cutting 3/8 Circle 30mm
SN764G	A664N	Nylon 2-0 Black 18" FS Reverse Cutting 3/8 Circle 26mm
SP1634	J8634N	Polypropylene 4-0 Blue 18" PC-3 Precision Conventional Cutting 3/8 Circle 16mm
SP1635	J8635N	Polypropylene 5-0 Blue 18" PC-3 Precision Conventional Cutting 3/8 Circle 16mm
SP1643	J8632N	Polypropylene 3-0 Blue 18" PC-5 Precision Conventional Cutting 3/8 Circle 19mm
SP1644	J8631N	Polypropylene 4-0 Blue 18" PC-5 Precision Conventional Cutting 3/8 Circle 19mm
SP5663G	J8663N	Polypropylene 3-0 Blue 18" PS-1 Precision Reverse Cutting 3/8 Circle 24mm
SP5680G	J8680N	Polypropylene 6-0 Blue 18" PS-3 Precision Reverse Cutting 3/8 Circle 16mm
SP5681G	J8681N	Polypropylene 5-0 Blue 18" PS-3 Precision Reverse Cutting 3/8 Circle 16mm
SP5682G	J8682N	Polypropylene 4-0 Blue 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SP5686G	J8686N	Polypropylene 5-0 Blue 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SP5687G	J8687N	Polypropylene 3-0 Blue 18" PS-2 Precision Reverse Cutting 3/8 Circle 19mm
SP5695G	J8695N	Polypropylene 6-0 Blue 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SP5697G	J8697N	Polypropylene 6-0 Blue 18" P-1 Precision Reverse Cutting 3/8 Circle 11mm
SP5698G	J8698N	Polypropylene 5-0 Blue 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SP5699	J8699N	Polypropylene 4-0 Blue 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SP5699G	J8699N	Polypropylene 4-0 Blue 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SP661G	J8661N	Polypropylene 5-0 Blue 18" FS-2 Reverse Cutting 3/8 Circle 19mm
SP665G	J8665N	Polypropylene 3-0 Blue 18" FS-2 Reverse Cutting 3/8 Circle 19mm
SP683G	J8683N	Polypropylene 4-0 Blue 18" FS-2 Reverse Cutting 3/8 Circle 19mm
SP684G	J8684N	Polypropylene 3-0 Blue 18" FS-1 Reverse Cutting 3/8 Circle 24mm
SP868G	J8685N	Polypropylene 2-0 Blue 18" FS Reverse Cutting 3/8 Circle 26mm
SS1639G	D639N	Silk 6-0 Black 18" P-1 Precision Reverse Cutting 3/8 Circle 11mm
SS5640	D640N	Silk 5-0 Black 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SS5641G	D641N	Silk 4-0 Black 18" P-3 Precision Reverse Cutting 3/8 Circle 13mm
SS5684	D1684N	Silk 3-0 Black 18" PS-1 Precision Reverse Cutting 3/8 Circle 24mm
SS678	D580N	Silk 0 Black 18" PSL Reverse Cutting 3/8 Circle 30mm

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
SS678	D678N	Silk 0 Black 18" FSL Reverse Cutting 3/8 Circle 30mm
SS684G	D684N	Silk 3-0 Black 18" FS-1 Reverse Cutting 3/8 Circle 24mm
SS685G	D685N	Silk 2-0 Black 18" FS Reverse Cutting 3/8 Circle 26mm
SS70M	D590N	Silk 0 Black 30" PSL Reverse Cutting 3/8 Circle 30mm
SS70M	D680N	Silk 0 Black 30" FSL Reverse Cutting 3/8 Circle 30mm
VP205X	J8554D	Polypropylene 5-0 Blue 24" RB-2 Taperpoint 1/2 Circle 13mm Double Armed
VP37MX	J8755N	Polypropylene 7-0 Blue 30" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP422X	J8522D	Polypropylene 3-0 Blue 36" SH Taperpoint 1/2 Circle 26mm Double Armed
VP521X	J8521D	Polypropylene 4-0 Blue 36" SH Taperpoint 1/2 Circle 26mm Double Armed
VP523	J8523D	Polypropylene 2-0 Blue 36" SH Taperpoint 1/2 Circle 26mm Double Armed
VP542X	J8534D	Polypropylene 3-0 Blue 48" SH Taperpoint 1/2 Circle 26mm Double Armed
VP543	J8533D	Polypropylene 2-0 Blue 48" SH Taperpoint 1/2 Circle 26mm Double Armed
VP556X	J8556D	Polypropylene 5-0 Blue 36" RB-1 Taperpoint 1/2 Circle 17mm Double Armed
VP557X	J8557D	Polypropylene 4-0 Blue 36" RB-1 Taperpoint 1/2 Circle 17mm Double Armed
VP558X	J8558D	Polypropylene 3-0 Blue 36" RB-1 Taperpoint 1/2 Circle 17mm Double Armed
VP566X	J8555D	Polypropylene 5-0 24" RB-1 RB-1 Taper Point 1/2 Circle 17mm
VP580X	J8580D	Polypropylene 5-0 Blue 36" BB Taperpoint 3/8 Circle 17mm Double Armed
VP581X	J8581D	Polypropylene 4-0 Blue 36" BB Taperpoint 3/8 Circle 17mm Double Armed
VP701X	J8701D	Polypropylene 7-0 Blue 18" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP702X	J8702D	Polypropylene 7-0 Blue 24" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP703X	J8805D	Polypropylene 6-0 Blue 24" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP706X	J8706D	Polypropylene 6-0 Blue 30" C-1 Taperpoint 3/8 Circle 13mm Double Armed
VP709X	J8709D	Polypropylene 6-0 Blue 30" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP710X	J8710D	Polypropylene 5-0 Blue 30" RB-2 Taperpoint 1/2 Circle 13mm Double Armed
VP713X	J8703D	Polypropylene 7-0 Blue 30" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP714X	J8714D	Polypropylene 6-0 18" RB-2 RB-2 Taper Point 1/2 Circle 13mm
VP715X	J8713D	Polypropylene 5-0 18" RB-2 RB-2 Taper Point 1/2 Circle 13mm
VP718X	J8718D	Polypropylene 6-0 Blue 18" C-1 Taperpoint 3/8 Circle 13mm Double Armed
VP719X	J8717D	Polypropylene 5-0 Blue 18" C-1 Taperpoint 3/8 Circle 13mm Double Armed
VP720X	J8720D	Polypropylene 5-0 Blue 36" C-1 Taperpoint 3/8 Circle 13mm Double Armed
VP725X	J8725D	Polypropylene 5-0 Blue 24" C-1 Taperpoint 3/8 Circle 13mm Double Armed
VP726X	J8726D	Polypropylene 6-0 Blue 24" C-1 Taperpoint 3/8 Circle 13mm Double Armed
VP72MX	J8747N	Polypropylene 7-0 Blue 24" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP733X	J8711D	Polypropylene 6-0 Blue 30" RB-2 Taperpoint 1/2 Circle 13mm Double Armed
VP746X	J8806D	Polypropylene 6-0 Blue 18" BV-1 Taperpoint 3/8 Circle 9.3mm Double Armed
VP74X	J8610N	Polypropylene 6-0 Blue 24" BV Taperpoint 3/8 Circle 11mm Double Armed

N=box of 12; D=box of 36

Covidien Cross Reference Guide

Covidien	SSC	Product Description
VP760X	J8762D	Polypropylene 3-0 Blue 30" SH-1 Taperpoint 1/2 Circle 22mm Double Armed
VP776X	J8776N	Polypropylene 6-0 Blue 30" BV Taperpoint 3/8 Circle 11mm Double Armed
VP832X	J8832N	Polypropylene 3-0 Blue 30" SH Taperpoint 1/2 Circle 26mm
VP833	J8833N	Polypropylene 2-0 Blue 30" SH Taperpoint 1/2 Circle 26mm
VP834	J8418N	Polypropylene 0 Blue 30" MO-6 Taperpoint 1/2 Circle 26mm
VP834	J8834N	Polypropylene 0 Blue 30" SH Taperpoint 1/2 Circle 26mm
VP842X	J8842D	Polypropylene 3-0 Blue 36" MH Taperpoint 1/2 Circle 36mm Double Armed
VP843	J8843D	Polypropylene 2-0 Blue 36" MH Taperpoint 1/2 Circle 36mm Double Armed
VP853	J8853D	Polypropylene 2-0 48" MH MH Taper Point 1/2 Circle 36mm
VPF543	J8533D	Polypropylene 2-0 Blue 48" SH Taperpoint 1/2 Circle 26mm Double Armed
VPF557X	J8557D	Polypropylene 4-0 Blue 36" RB-1 Taperpoint 1/2 Circle 17mm Double Armed

N=box of 12; D=box of 36



Order Online www.SurgicalSpecialties.com



Order by Phone **877-991-1110**

Sharpoint™ is also available through major healthcare distributors.

**SURGICAL
SPECIALTIES**
CORPORATION™

For More Information, Visit SurgicalSpecialties.com

Surgical Specialties Corporation, 247 Station Drive, Suite NE1, Westwood, MA 02090, U.S.A.
Phone: 610.404.3520 • 877.991.1110 (USA and Canada only) • Fax: 610.404.4010
All trademarks are property of their respective owners. ©2018 Surgical Specialties Corporation.
All Rights Reserved. SSC-028 R5 03/2020