

## Simplify skin closure with Quill® barbed suture.



### What is Quill barbed suture?

It's a proven wound closure device that eliminates the need to tie knots during soft tissue approximation.

Quill uses a running stitch to hold wound edges together with the strength of an interrupted suture.

The barbs lock the suture in place, maintaining tension the length of the wound without assistance.

Variable loop provides  
a simple fixation method.



Retractable barbs  
allow the suture to  
pass smoothly through  
tissue in one direction.

Barbs anchor  
into the tissue for  
secure fixation.

### WHY QUILL?

**Produces better apposition.**

**Improves cosmesis.**

**Eliminates knots for a faster close.**

**Requires less foreign material,  
no knots under the skin.**

**Avoids irritating staples.**

### Free sampling program for evaluation

See how quickly and effectively you can  
close wounds with no knots.

**Learn from Corza Medical's suture specialist —**  
virtual and hands-on training and support.

**Need to test the waters?** We will give you  
**FREE TruStitch™** to stitch over Quill while  
you learn to go knotless.



## SELECT FROM OUR TOP 20 SKUs:

Monoderm™				USP Size (Box of 12)			
Reverse Cutting Needles ▽		Color	Length	3-0	2-0	0	2
FS-2	19mm 3/8 Circle	Violet	30cm 60cm	YA033R5 YA036R5			
FS-1	24mm 3/8 Circle	Violet	30cm 60cm	YA033R6	YA023R6 YA026R6		
FS	26mm 3/8 Circle	Violet	60cm		YA026R7		
FSL	30mm 3/8 Circle	Violet	30cm			YA03R8	
CP-2	26mm 1/2 Circle	Violet	30cm			YA03R14	
Taper Point Needles ●		Color	Length	3-0	2-0	0	
RB-1	17mm 1/2 Circle	Violet	30cm		YA023T2		
SH	26mm 1/2 Circle	Violet	30cm		YA023T4		

PDO (Polydioxanone)				USP Size (Box of 12)			
Reverse Cutting Needles ▽		Color	Length	3-0	2-0	0	2
FS-I	24mm 3/8 Circle	Violet	30cm	RA033R6	RA023R6		
FSL	30mm 3/8 Circle	Violet	30cm			RA03R8	
CP-2	26mm 1/2 Circle	Violet	30cm			RA03R14	
CP	40mm 1/2 Circle	Violet	60cm				RA26R16
Bi-Directional Reverse Cutting Needles ▽		Color	Length	3-0	2-0	0	2
FS-I	24mm 1/2 Circle		14x14	RDA031R6	RDA021R6		
Taper Point Needles ⊙		Color	Length	3-0	2-0	0	2
CT-2	26mm 1/2 Circle	Violet	30cm			RA03T6	
CT-I	36mm 1/2 Circle	Violet	30cm 60cm				RA23T7 RA26T7

View our full catalog at [www.surgicalspecialties.com/quill-vet-kit/](http://www.surgicalspecialties.com/quill-vet-kit/)

### SIMPLIFIED PRICING

Standardized pricing across materials regardless of needle type, size, length, color.

MATERIAL TYPE	QUILL LIST PRICE UNIDIRECTIONAL	QUILL LIST PRICE BIDIRECTIONAL
MONODERM™	\$125.00	\$145.00
PDO™	\$136.00	\$158.00

**CONTACT YOUR CORZA  
SUTURE SPECIALIST TODAY:**

**Email** [animalhealth@corza.com](mailto:animalhealth@corza.com)  
to select your **FREE** samples  
and kick off your evaluation!